

Analytical Balance - Auto-CAL

ATN-G Series | HSN Code: 90273020



animals weighing



density



AutoCAL



recipe



statistics



pieces counting



percentage



summation



date/time



threshold/selection



max weight



min weight



products and users database

Option:-



alibi memory



measurements database

FEATURES

- Complete automatic time & temperature calibration with GMP-GLP record printout
- Brilliant Graphic LCD display with simple navigation windows and menu to guide you through balance functions
- User selectable filter, stability indicator and auto-zero tracking for use in difficult environment
- Procell software allows for direct transfer of weighing result to excel spreadsheet
- Weigh below hook for below weighing application with in-built density measurement software
- With a die-cast metal base and stainless-steel pan, and balance is built to withstand rugged every day Weighing
- RS232C and USB connector to connect with PC

Application:

- Parts Counting
- Animal Weighing
- Check-Weighing
- Display Hold
- Percent Weighing
- Totalizing
- SQC Quality Control

MODEL	CAPACITY (gm)	RESOLUTION (gm)	PAN mm	REPEATABILITY (gm)	LINEARITY (gm)	RESPONSE TIME (Sec.)
Auto-Internal Calibration						
ATN 120G	120	0.0001	Ø 90	± 0.0001	± 0.0002	< 3
ATN 160G	160	0.0001	Ø 90	± 0.0001	± 0.0002	< 3
ATN 220G	220	0.0001	Ø 90	± 0.0001	± 0.0002	< 3

Optional

• Density Kit



• Data Printer



• Software



• Test Weights



Product Range

• pH/mV/°C/ORP Analyzer • Conductivity/TDS/Resistivity/Salinity/°C Analyzer
 • Multiparameter Analyzer (pH/ORP/Conductivity/TDS/Resistivity/Salinity/°C)
 • Karl Fischer Titrator • Auto Titrator • Automatic Melting Point Apparatus
 • Flame Photometer • Friability Test Apparatus • Tablet Dis-integration Test Apparatus
 • Tap Density Test Apparatus • Tablet Hardness Tester • Leak Test Apparatus
 • Tablet Dissolution Test Apparatus • Electromagnetic Sieve Shaker • Ultrasonic Bath (Sonicator)
 • Analytical Electronics Balance • Halogen Based Moisture Analyzer
 • Visible & UV-Vis. Spectrophotometer • Nitrogen Gas Generator • Zero Air Gas Generator
 • Combination Nitrogen / Zero Air Gas Generator • Air Drying Unit for TOC Analyzer
 • Dew Point Apparatus • Gas Handling & Purification System • Laboratory Furniture & Fume Hood
 • Exhaust Hood for AAS/ICP.

ANALAB SCIENTIFIC INSTRUMENTS PVT. LTD.

AN ISO 9001 : 2015 CERTIFIED CO.
 CIN: U33110GJ2005PTC046555

Plot No. 361, Sayajimarket,
 N. H. No. 8 By-Pass,
 At & Post : Sayajipura,
 Vadodara - 390019, Gujarat, India.

+91 265 2514404, 9714503393

info@analab.co.in

www.analab.co.in, www.analabscientific.com



Specification & design are subject to change without prior notice for our continuous development.